AT04509

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The AN TAR Chest-Rib Belt for Men (AT04509) is made of soft-touch elastic material for a perfect fit and compression. Thanks to the hook-and-loop fastening and the side supports, the product has a wide adjustment range and provides adequate stabilisation.

TARGET PATIENT GROUP

Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see section on intended use of the device in these instructions). The device can be purchased by the user on their own or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions/size/variants of the device, the indications and contraindications for use, as well as the information provided by the manufacturer.

CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) that prevent safe handling of the product.

INDICATIONS/APPLICATIONS:

- after rib and sternum fractures
- after cardiac surgery
- after thoracic surgery

FEATURES:

- Reduces pain and stabilises the chest
- Improves post-operative wound healing conditions
- Prevents post-operative wound dehiscence
- Provides comfort and facilitates recovery

APPLICATION: Unfasten the Velcro fasteners. Then place the belt on the chest by fastening the Velcro at the front. Adjust the length of the braces with the Velcro.

SIZE - Chest circumference (cm)

S (75-84cm)
M (85-94)
L (95-104)
XL (105-114)
XXL (115-125)
XXXL (126-135)

WASHING INSTRUCTIONS/MAINTENANCE

Fasten the Velcro fasteners before washing. Wash by hand in water up to 40 °C with gentle detergents, do not use strong detergents. Do not iron or bleach. Dry flat in the shade, avoid direct contact with the sun or other heat source. Do not wring out - squeeze out excess water by towel drying.

NOTES:

- Read the contents of this manual before wearing the belt.
- Before each use, make sure there is no abrasion of the seams and individual components on the belt.
- Use of the product contrary to the instructions may result in serious complications.
- The belt should be used as recommended by a doctor or physiotherapist.
- The product should not be used in the case of open wounds, skin abrasions or after the application of warming gel or ointment.
- Wearing the strap overnight is only possible if recommended by the attending physician.
- Unless otherwise advised by the doctor, the belt should not be worn for more than 4 hours a day.
- If acute pain or skin irritation occurs while wearing the product, remove the product immediately and consult your doctor.

STORAGE

Store the device in a dry cool room, away from moisture and direct sunlight.

CAUTION:

In the event of a device-related „serious incident“ that has directly or indirectly led to, could have led to or is likely to lead to any of the following:

- a) death of a patient, user or other person, or
- b) temporary or permanent deterioration of the health of a patient, user or other person, or
- c) a serious risk to public health

the above ‚serious incident‘ must be notified to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTES:

In the event of pain, allergic reactions or other disturbing symptoms unclear to the user related to the use of the medical device, consult a healthcare professional.

HOW TO DISPOSE OF THE DEVICE AT END-OF-LIFE

Once the device has been withdrawn from service, the medical device can be disposed of as normal municipal waste.

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

 ANTAR Sp. J.
03-068 Warsaw, ul. Zawiańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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