AT53060, AT53061
ORTEZA KOLANA INSTRUKCJA OBSŁUGI
KNEE ORTHOSIS INSTRUCTIONS FOR USE
GEBRAUCHSANWEISUNG FÜR DIE KNIEORTHESE
NÁVOD K POUŽITÍ KOLENNÍ ORTÉZY
NÁVOD NA POUŽITIE - ORTÉZA NA KOLENO
ORTHÈSE DU GENOU MODE D'EMPLOI
GEBRUIKSAANWIJZING KNIEORTHESE
ORTESIS DE RODILLA INSTRUCCIONES DE USO
ISTRUZIONI PER L'USO DELL'ORTESI DI GINOCCHIO
KNÄORTOS BRUKSANVISNING

MD

PL

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ENG

GENERAL WARNING

Thank you for your purchase. Please read and understand these instructions before using the product. If any warnings or instructions are unclear, contact your healthcare professional or dealer to prevent product damage.

ATTENTION: It is forbidden to use the product other than for its intended purpose.

ATTENTION: Before use, check that no part is damaged and looks correct.

ATTENTION: The manufacturer isn't liable for damage due to negligent maintenance, insufficient servicing, or failure to follow the manual's instructions.

CAUTION:

In case of a device-related "serious incident" leading to or likely to lead to:

- (a) death of a patient, user or other person, or
- (b) temporary or permanent deterioration of the health of a patient, user or other person, or
- (c) a serious risk to public health

Report the mentioned 'serious incident' to both the manufacturer and the relevant authority in the user or patient's residing Member State. In Poland, contact the Office for Registration of Medicinal Products, Medical Devices, and Biocidal Products.

NOTES:

If experiencing pain, allergies, or any unclear distress while using the medical device, consult a healthcare professional.

TARGET PATIENT GROUP

People with conditions that align with the device's intended use (as outlined in the instructions) can buy it independently or through a doctor's, therapist's, or specialist's recommendation. Consider the device's available sizes, necessary functions, variants, indications, contraindications, and the manufacturer's provided information, regardless of how it's purchased.

SCOPE OF SUPPLY

Orthosis, instructions for use

USE/ INDICATIONS

- acute/chronic pain syndromes
- treatment of knee joint injuries requiring limitation of the range of motion
- injuries and dysfunctions without the need for immobilisation of the joint
- support of the ligamentous system of the knee joint

CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) which prevent safe handling of the product.

DISPOSAL OF THE PRODUCT

The product must be disposed of in accordance with the applicable disposal conditions.

CLEANING/MAINTENANCE

Hand wash the product at temperatures up to 40°C using a mild detergent, then allow it to air dry. Avoid using harsh objects, solvents, strong detergents, tumble drying, dry cleaning, or bleach.

STORAGE

Store the product in a dry cool room, away from moisture and direct sunlight.

FUNCTIONS AND FEATURES OF THE PRODUCT

- Ability to stabilise and lock the knee joint at a specific angle
- Possibility to adjust the range of flexion of the knee joint
- Wide range of movement
- AT53061 - Supports made of carbon fibre

SIZES

AT53060 - UNI size: 49 cm long

AT53061 - UNI size: length adjustable 46-66 cm

PARAMETERS AT53060

Weight: 0.74 kg

Flexion settings: 0, 10, 20, 30, 45, 60, 75, 90, 105, 120°

Extension settings: 0, 10, 20, 30°

PARAMETERS AT53060

Weight - 0.64 kg

Flexion settings: -10, 0, 10, 20, 30, 45, 60, 75, 90, 105, 120°

Upright settings: -10, 0, 10, 20, 30°

INSTRUCTIONS FOR USE AT53060

NOTE: The product should be fitted and adjusted with a doctor or specialist present on a slightly bent or straight leg.

1. Unfasten and remove the Velcro straps from the plastic buckles and then loosen the straps (Fig. 1a).
2. put the brace on the leg, the side rails should be on both sides of the leg and the clocks should be in the middle of the knee (Fig. 2a).
3. Secure the orthosis by holding it with your hand to prevent movement while fastening. Attach the buckles on the handles and adjust the stabilization level by altering the strap length. Fasten the Velcro straps in sequence, following the order shown in Figures 4a to 4c.

Fasten the straps in the following order:

Strap 1: Fasten the strap on the lower part of the thigh.

Strap 2: Fasten the strap on the upper part of the calf.

Strap 3: Fasten the strap on the upper thigh.

Strap 4: Fasten the strap on the lower part of the calf. Fastening too tightly can lead to slowed blood circulation or a downward displacement of the bracket. Fastening too loosely will cause the orthosis to fall off.

Note: Position the hinge slightly behind the knee's midline when viewed from the side. Modify the strap positions if needed. After adjusting, the orthosis can be removed using the buckles without altering each attached strap (refer to Figure 4a).

The orthosis is adaptable for use with a locked knee position or with flexion. Choose the desired flexion angle on the clock to set the range of adjustment. Set the correct flexion angle as shown in Figure 5. Additionally, you can lock the orthosis by fixing the desired flexion angle.

INSTRUCTIONS FOR USE AT53061

NOTE: The product should be fitted and adjusted in the presence of a doctor or specialist present on a slightly bent or straight leg.

1. Loosen the straps and release the buckles one at a time from the handles, then loosen the straps further by sliding the buckle (Fig. 1a).
2. put the orthosis on your leg, the side rails should be on both sides of the leg and the clocks should be at the mid-knee level (Fig. 2a).
3. adjust the orthosis length by pressing the telescopic adjustment button on the rail, moving it to the desired length, and then releasing the button. Ensure the splint is securely locked by placing all four buttons in the correct locking position during fitting.

IMPORTANT:

Doctors or specialists should adjust the orthosis length for optimal knee stability. Position the splint height to end 5 cm below the groin and ensure the lower part of the splint doesn't reach the Achilles tendon. NOTE: The product should be worn on a slightly bent or straight leg.

Hold the brace to prevent movement and pass the Velcro-ended straps through the plastic buckles. Adjust the Velcro straps to fit the brace around your leg, fastening the straps in sequence as indicated (Fig. 4a-4c).

Fasten the straps in the following order:

Strap 1: Fasten the strap on the lower part of the thigh.


Strap 2: Fasten the strap on the upper part of the calf.

Strap 3: Fasten the strap on the upper thigh.

Strap 4: Fasten the strap on the lower part of the calf. Secure the strap on the upper calf but avoid fastening it too tightly to prevent restricted blood flow or splint displacement. Conversely, if it's fastened too loosely, the orthosis might fall off.

Note: Position the hinge slightly behind the knee's midline when viewed from the side. Adjust strap positions as needed. After adjusting, the orthosis can be removed using the buckles without the need to loosen or tighten each attached strap (refer to Figure 4a).

The orthosis can function in a locked knee position or with flexion. To adjust the flexion angle, unlock the fixed position by sliding the lock button on the hinge and set the desired flexion angle (refer to Figure 5). By moving the lock button on the hinge, the orthosis can be secured at the desired angle.

IMPORTANT:

- 1) The products are designed for individual patients and are not recommended for multiple patients.
- 2) Both hinges must be set to the same upright and/or flexion setting.
- 3) The patient should consult a doctor before adjusting the orthosis.
- 4) Proper rehabilitation and activity improvement are also essential processes in safe treatment planning.

STRAP LENGTH ADJUSTMENT AT53061

Measure the needed length, open the Velcro at the strap end, detach the Velcro, trim the strap to the desired length, and reattach the Velcro at the strap end using glue.

LOCKING INSTRUCTIONS AT53061

- 1) Determine the desired locking angle.
- 2) Move the adjustment button towards the centre of the hinge to lock the orthosis in the correct hinge angle position.
- 3) To adjust within the 0-30° range, slide the locking button outward. Pull the EXTENSION button to reach the 0° position and the FLEXION button to attain the 30° position.
- 4) Repeat the process to set a different range, like 30-120°. Pull the EXTENSION button to 30° and the FLEXION button to 120°.

LOCKING INSTRUCTIONS AT53060

- 1) Specify the desired locking angle.
- 2) Adjust to a 0-30° range by pulling the EXTENSION button to 0° and the FLEXION button to 30°.
- 3) To configure a different range, such as 30-120°, adjust by pulling the EXTENSION button to 30° and the FLEXION button to 120°.
- 4) To lock the orthosis in the chosen position (30/30), position the EXTENSION button at 30° and the FLEXION button at 30°.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. AN TAR Sp. J., 03-068 Warsaw, ul. Zawiślanska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period AN TAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. AN TAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to AN TAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślanska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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