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ORTEZA KOLANA: AT53010, AT53011, AT53012, AT53013.

ORTEZA STAWU SKOKOWEGO: AT53014, AT53015, AT53016, AT53017, AT92204

ORTEZA NADGARSTKA: AT53018, AT53019, AT53020, AT53021, AT53055 ORTEZA ŁOKCIA: AT53022, AT53023, AT53024, AT53025

INSTRUCTIONS FOR USE

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GENERAL WARNING

Do not use the product before reading and understanding these instructions. If you do not understand the warnings, cautions or instructions, contact your healthcare professional or dealer to avoid damage to the product.

APPLICATION/ INDICATIONS

Wrist joint orthoses:

- prevention of injuries and strains to the wrist joint
- protection of the injured joint
- stabilise the joint, reduce the risk of aggravation of the injury

Elbow joint orthoses:

- prevention of injuries and overloading of the elbow joint
- protection of the injured joint
- stabilisation of the joint, reducing the risk of injury aggravation

Knee joint orthoses:

- prevention of injuries and overloading of the knee joint
- protection of the injured joint
- stabilise the joint, reduce the risk of aggravation of the injury

Ankle orthoses:

- prevention of injuries and overloading of the ankle joint
- protection of the injured joint
- stabilise the joint, reduce the risk of aggravating the injury
- support for the treatment of sprains and dislocations of the ankle joint

CONTRAINDICATIONS

Physical or mental limitations (e.g. visual impairment) that prevent safe handling. Open wounds, abrasions of the epidermis, and previous use of a warming gel or ointment.

NOTES:

- Read the contents of this manual before putting on the orthosis.
- Before each use, make sure there is no abrasion of the seams and individual components on the orthosis.
- Use of the product contrary to the instructions may result in serious complications.
- Use the brace as recommended by your doctor or physiotherapist.
- The orthosis may only be worn overnight if recommended by the treating doctor.

SIZES

Wrist orthoses: to determine the size, measure the circumference around the wrist

AT53021, AT53019 - (UNI)

AT53018, AT53020 (S 16-20 cm, M 21-24 cm, L 25-29 cm, XL 30-34 cm)

AT53055 (S 11-14 cm, M 14-17 cm, L 17-20 cm, XL 20-25 cm)

Elbow orthoses: to determine the size, measure the circumference around the elbow

AT53022, AT53024, AT53025 (S 20-23 cm, M 23-26 cm, L 26-29 cm, XL 29-32 cm)

AT53023 - (UNI)

Knee bracing: to determine the size, measure the circumference around the knee at 20° of flexion

AT53011 - (UNI)

AT53012, AT53013, AT53010 (S 29-34 cm, M 35-39 cm, L 40-49 cm, XL 45-49 cm)

Ankle orthosis: (metatarsal circumference:) AT53015 - (UNI)

AT53014, AT53016, AT53017 (S 10-18 cm, M 18,5-23 cm, L 23-33 cm, XL:28-38 cm)

AT92204 - (S 15.6 cm - 21.0 cm, M 21.6 cm - 26.0 cm, L 26.7 cm - 31.1 cm, XL 31.8 cm - 36.3 cm,

XXL 36.8 cm - 41.3 cm)

CLEANING/MAINTENANCE:

- The brace can only be washed by hand in lukewarm water (up to 40°C) with a mild detergent (e.g. grey soap).
- Do not iron.
- Do not dry-clean.
- Do not use bleach.
- Protect from moisture.
- Do not expose to high temperatures.
- Do not tumble dry.
- Dry in the shade when unfolded.

STORAGE

Store the product in a dry cool room, away from moisture and direct sunlight.

CAUTION:

In the event of a device-related "serious incident" that has directly or indirectly led to, could have led to or is likely to lead to any of the following:

- (a) death of a patient, user or other person, or
- (b) temporary or permanent deterioration of the health of a patient, user or other person, or
- (c) a serious risk to public health

the above 'serious incident' must be notified to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. In the case of Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

NOTES:

In case of pain, allergic reactions or other disturbing symptoms unclear to the user related to the use of the medical device, consult a healthcare professional.

CAUTION:

Use of the product other than for its intended purpose is prohibited.

CAUTION:

The manufacturer is not responsible for damage caused by negligence of maintenance, cleaning rules or resulting from non-compliance with the recommendations contained in this manual

SCOPE OF SUPPLY

Orthosis, instructions for use.

TARGET PATIENT GROUP

Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see section on intended use of the device in this instruction manual). The device can be purchased by the user on their own or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions and variants of the device, the indications and contraindications for use and the information provided by the manufacturer.

DISPOSAL OF THE PRODUCT

The product must be disposed of in accordance with the applicable disposal conditions.

AT53012, AT53016, AT53020, AT53024: 15% nylon, 2% spandex, 25% viscose, 33% cotton, 25% latex thread Note: viscose is made from bamboo

AT53011, AT53015, AT53019, AT53023: 65% polyester, 35% elastic thread

AT53010, AT53014, AT53018, AT53022: 70% polyester, 5% spandex, 25% latex thread

AT53013, AT53017, AT53021, AT53025, AT53055: 70% nylon, 5% spandex, 25% latex thread

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WARRANTY CARD

Model:
Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR MEDICAL Sp. z o.o., 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR MEDICAL Sp. z o.o.
03-068 Warsaw, ul. Zawisłańska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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