AT51424

**INSTRUKCJA UŻYTKOWANIA, KOPUŁA MASUJĄCO-REHABILITACYJNA
USER MANUAL, MASSAGE AND REHABILITATION DOME
GEBRAUCHSANWEISUNG, MASSAGE- UND REHABILITATIONSKUPPEL
NÁVOD K POUŽITÍ, MASÁŽNÍ A REHABILITAČNÍ KOPULE
NÁVOD NA POUŽITIE, MASÁŽNA A REHABILITAČNÁ KOPULA
MODE D'EMPLOI, DÔME DE MASSAGE ET DE RÉÉDUCATION
GEBRUIKSAANWIJZING, MASSAGE- EN REVALIDATIEDOME
INSTRUCCIONES DE USO, CÚPULA DE MASAJE Y REHABILITACIÓN
ISTRUZIONI PER L'USO, CUPOLA DA MASSAGGIO E RIABILITAZIONE
BRUKSANVISNING, MASSAGE- OCH REHABILITERINGSKUPOL**

MD

PL

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EN

Thank you for purchasing our product. Please read the instruction manual carefully. Do not use the product before reading and understanding these instructions. If you do not understand the warnings, notes, or recommendations, please consult a healthcare specialist or the seller to avoid product damage or bodily injury.

WARNING: Check all parts for any damage incurred during transportation. If such damage is observed, DO NOT use the product. More information is available from the seller.

PURPOSE/INDICATIONS

The balance ball is a multifunctional tool widely used in rehabilitation. It is designed to improve ankle stability, enhance shoulder girdle endurance, and provide sensory receptor massage for feet and hands. With its unique properties, it helps improve body stability and coordination, making it ideal for balance training, including fitness activities and floor exercises. Specially designed spikes stimulate circulation and provide sensory receptor massage, adding to the comfort of use. For better training results, purchasing multiple balance balls is recommended, enabling the creation of more advanced exercise circuits.

Product Features:

- The set includes 2 air-filled balance balls.
- Supports ankle joint stabilization.
- Enhances shoulder girdle muscle endurance.
- Improves motor coordination and proprioception.
- Textured, spiked surface provides optimal sensory stimulation.
- Easy to use with adjustable pressure levels.
- Suitable for all age groups and skill levels.
- Made without latex or phthalates, safe for allergic individuals.

Technical Specifications:

- Dome diameter: 16 cm
- Color: blue
- Maximum user weight: 180 kg
- Composition: PVC

This symbol indicates the maximum user weight.

CONTRAINDICATIONS

Physical or mental limitations (e.g., visual impairment) that prevent safe handling of the product.

CONTENTS OF DELIVERY

Massage and rehabilitation dome, instruction manual.

STORAGE

Store in a dry and cool place, away from direct sunlight.

MAINTENANCE

Clean the cushion with a damp cloth and mild detergent.

WARNING

Do not use strong detergents, sharp objects, or solvents for cleaning.

WARNING

In the event of a "serious incident" associated with the product that has directly or indirectly caused, could have caused, or may cause any of the following events:

- a) death of the patient, user, or another person,
- b) temporary or permanent deterioration in the health of the patient, user, or another person, or
- c) a serious threat to public health,

the above "serious incident" should be reported to the manufacturer and the competent authority of the member state where the user or patient resides. In Poland, the relevant authority is the Office for Registration of Medicinal Products, Medical Devices, and Biocidal Products.

WARNING

In case of pain, allergic reactions, or other concerning symptoms related to the use of the medical device, consult a healthcare professional.

WARNING

Using the product in any manner other than intended is prohibited.

WARNING

The manufacturer is not responsible for any damage caused by neglect of maintenance, improper servicing, or failure to follow the recommendations contained in this instruction manual.

WARNING: Flammable product.

WARNING: Product not suitable for children under the age of 3.

PRODUCT DISPOSAL AFTER END OF USE

Once out of use, the medical device can be disposed of as regular household waste.

TARGET PATIENT GROUP

Individuals dealing with diseases, dysfunctions, or injuries for the treatment/rehabilitation or compensation of which this product is intended (see the purpose section of this instruction manual). The product may be purchased independently by the user or based on recommendations from a doctor, therapist, or other specialist. Both for self-purchase and when recommended by a doctor/therapist/other specialist, it is essential to consider the available sizes/required functions and product variations, indications, and contraindications for use, as well as information provided by the manufacturer.

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Model:
Serial number:

Point of sale stamp

and a legible signature of the Seller:

1. AN TAR Sp. J., 03-068 Warsaw, ul. Zawiślńska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period AN TAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. AN TAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to AN TAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warsaw, ul. Zawiślńska 43
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net

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